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

Rycina 7

Bibliografia

- Goligher EC, Dres M, Fan E, Rubenfeld GD, et al. Mechanical Ventilation-induced Diaphragm Atrophy Strongly Impacts Clinical Outcomes. *Am J Respir Crit Care Med.* 2018 Jan 15;197(2):204-213.
- Dres M, Dubé BP, Mayaux J, et al. Coexistence and Impact of Limb Muscle and Diaphragm Weakness at Time of Liberation from Mechanical Ventilation in Medical Intensive Care Unit Patients. *Am J Respir Crit Care Med.* 2017 Jan 1;195(1):57-66.
- Dres M, Goligher EC, Heunks LMA, et al. Critical illness-associated diaphragm weakness. *Intensive Care Med.* 2017 Oct;43(10):1441-1452.
- Kim WY, Suh HJ, Hong SB, Koh Y, et al. Diaphragm dysfunction assessed by ultrasonography: influence on weaning from mechanical ventilation. *Crit Care Med.* 2011 Dec;39(12):2627-30.
- Colombo D, Cammarota G, Alemani M, et al. Efficacy of ventilator waveforms observation in detecting patient-ventilator asynchrony. *Crit Care Med.* 2011 Nov;39(11):2452-7.
- Goligher E. Diaphragm dysfunction: monitoring and mitigation during mechanical ventilation. Lecture recording, at 23:50: www.criticalcarenews.com/webinars-symposia/diaphragm-weakness-clinical-outcomes-mechanical-ventilation/
- Goligher EC, Schepens T. Using ultrasound to prevent diaphragm dysfunction. *ICU Management & Practice*, Volume 18 - Issue 4, 2018.
- Emeriaud G, Larouche A, Ducharme-Crevier L, et al. Evolution of inspiratory diaphragm activity in children over the course of the PICU stay. *Intensive Care Med.* 2014 Nov;40(11):1718-26.
- Bellani G, Pesenti A. Assessing effort and work of breathing. *Curr Opin Crit Care.* 2014 Jun;20(3):352-8.
- Bellani G, Coppadoro A, Patroniti N, et al. Clinical assessment of autospontaneous end-expiratory pressure by diaphragmatic electrical activity during pressure support and neurally adjusted ventilatory assist. *Anesthesiology.* 2014 Sep;121(3):563-71.
- Doorduyn J, Sinderby CA, Beck J, et al. Automated patient-ventilator interaction analysis during neurally adjusted noninvasive ventilation and pressure support ventilation in chronic obstructive pulmonary disease. *Crit Care.* 2014 Oct 13;18(5):550.
- Ducharme-Crevier L, Beck J, Essouri S, et al. Neurally adjusted ventilatory assist (NAVA) allows patient-ventilator synchrony during pediatric noninvasive ventilation: a cross-over physiological study. *Crit Care.* 2015 Feb 17;19:44.
- A Skorko, D Hadfi eld, A Vercueil, et al. Retrospective review of utilisation and outcomes of diaphragmatic EMG monitoring and neurally adjusted ventilatory assist in a central London teaching hospital over a 3-year period. *Critical Care* 2013, 17(Suppl 2):P146.
- Rahmani A, Ur Rehman N, Chedid F. Neurally adjusted ventilatory assist (NAVA) mode as an adjunct diagnostic tool in congenital central hypoventilation syndrome. *J Coll Physicians Surg Pak.* 2013 Feb;23(2):154-6.
- Stein H, Firestone K. Application of neurally adjusted ventilatory assist in neonates. *Semin Fetal Neonatal Med.* 2014 Feb;19(1):60-9.
- Thille AW, Rodriguez P, Cabello B, et al. Patient-ventilator asynchrony during assisted mechanical ventilation. *Intensive Care Med.* 2006 Oct;32(10):1515-22.
- Colombo D, Cammarota G, Bergamaschi V, et al. "Physiologic response to varying levels of pressure support and neurally adjusted ventilatory assist in patients with acute respiratory failure," *Intensive Care Medicine*, vol. 34, no. 11, pp. 2010-2018, 2008.
- Bellani G, Mauri T, Coppadoro A, et al. Estimation of patient's inspiratory effort from the electrical activity of the diaphragm. *Crit Care Med.* 2013 Jun;41(6):1483-91.
- Liu L, Liu H, Yang Y, et al. Neuroventilatory efficiency and extubation readiness in critically ill patients. *Crit Care* 2012; 16:R143.
- Jansen D, Jonkman AH, Roesthuis L, et al. Estimation of the diaphragm neuromuscular efficiency index in mechanically ventilated critically ill patients. *Crit Care.* 2018 Sep 27;22(1):238.
- Passath C, Takala J, Tuchscherer D, et al. Physiologic response to changing positive end-expiratory pressure during neurally adjusted ventilatory assist in sedated, critically ill adults. *Chest.* 2010 Sep;138(3):578-87.
- Barwing J, Pedroni C, Olgemöller U, et al. Electrical activity of the diaphragm (EAdi) as a monitoring parameter in difficult weaning from respirator: a pilot study. *Crit Care.* 2013 Aug 28;17(4):R182.
- Liu L, Liu H, Yang Y, et al. Neuroventilatory efficiency and extubation readiness in critically ill patients. *Crit Care*, 16 (2012), pp. R143.
- Rozé H, Repusseau B, Perrier V, et al. Neuro-ventilatory efficiency during weaning from mechanical ventilation using neurally adjusted ventilatory assist. *Br J Anaesth*, 111 (2013), pp. 955-960.

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Producent: Maquet Critical Care AB · 171 54 Solna, Sweden · Phone: +46 (0)10-335 00 00 · info@getinge.com